


Detergent for automated washing process of animal cages Liquid concentrate

Fields of application

Automated cleaning of animal rodent cages, bottles and accessories as well as cage, storage and transport racks in specially designed washers.



Available in 10 and 20 L

Performance Spectrum

- Ecologically friendly neutral detergent.
- Effective for the removal of light organic soiling associated with small rodent caging.
- Compatible with Stainless Steel, Polycarbonate, Polyetherimide, Polysulfone and Polyphenylsulfone.

Special Properties

- Free of phosphates.
- Without oxidizing agents.
- Excellent material compatibility.

Application and dosage

Allentown Neutral detergent is ecologically friendly for use within industrial cage washers for the effective washing, mice and rat caging systems and associated components. Dose rates are determined by the level of soiling, cycle configuration and water conditions, and is recommended for use between 0.35-0.5 % (3.5-5 mL/L).

Compatible for all types of polysulfone and polycarbonate caging, enrichment and stainless steel components. Suitable for use with automatic dosing systems and bulk delivery systems.

Optimal water conditions for use with the detergent should be between $\text{CaCO}_3 > 50-100$ ppm.

General instructions on use

- Allentown Neutral detergent does not require neutralisation.
- Do not mix with other products.
- Purge dosing system including suction hose with water before changing product.
- The instructions given by the manufacturer of the washers are to be observed.
- Please observe the recommendations on cage washing given by the manufacturer of the items to be washed and the recommendations of the working group for cage processing (AK KAB) in the current issue of the AK KAB brochure "Cage Processing in Animal Facilities properly done".
- For professional use only.

Dosage Guide		
Rodent Caging (Mice & Rats)	0.35-0.5 %	3.5-5 mL per litre of water

Technical data

pH-range: 7

Ingredients

Ingredients are compliant to the provisions of Regulation (EC) No 648/2004 on detergents:
10-20 % Sodium Hydrogencarbonate.

Storage information

Always store at a temperature ranging between 0 and 30 °C.

Usable for 4 years from manufacture date, when stored as recommended. Please see container for dedicated batch details.

Hazard warnings and safety advice

The concentrate is classified and labelled according to the CLP Regulation (EC) No 1272/2008.

This mixture is non-hazardous.

Precautionary Statements:

- Do not ingest product.
- It is recommend to practice good hygiene measures when using this product, wash hands after use.

For further safety information see Safety data sheets.