


# RoVent® Jr.

Small animal automated ventilator



## Features & Benefits:

- **Free Offer:**  
RightTemp® temperature monitor and homeothermic control module  
*Warm and auto-regulate the animal's body temperature*
- **Fully automatic with touchscreen control**  
*Safely ventilate using auto-calculated settings based on animal weight*
- **Ventilate animals from 3g to 1,250g**  
*Ideal for mice, rats and small animals*
- **Volume control and pressure control**  
*Reduce footprint with one, dual control ventilator*
- **Multi-animal capability\***  
*Greater efficiency*

*\*Our clinical department will be pleased to review your study protocol to confirm that our product will meet your specific study requirements.*

1-Year 100% Parts & Labor Warranty (3-Year extended warranty available)

The RoVent Jr. small animal automated ventilator is ideal for mice, rats and small animals up to 1,250 grams.

RoVent Jr. is fully automatic — simply enter the animal's weight and press run.

[kentscientific.com](http://kentscientific.com)

Order your RoVent Jr.:  
[www.animalab.eu](http://www.animalab.eu) | [info@animalab.eu](mailto:info@animalab.eu)

Choose the ventilator that fits your research needs:

## Comparison Chart

### RoVent® Jr. & PhysioSuite® with RoVent Module



	RoVent Jr.	PhysioSuite with RoVent Module
Volume ventilation	●	●
Pressure ventilation	●	●
Pulse ventilation mode	○	●
PEEP	○	●
Assist mode	○	●
Sigh breath	○	●
I:E ratio	●	●
Tidal volume – measured	●	●
Auto set-up	Limited	●
User profiles	○	●
Alarms (parameter)	○	●
Trigger in	○	●
Analog out	○	●
Digital out	○	●
Sync out	○	●
ADI integration	○	●
Data output	○	●
Internal data logging	○	●
Realtime clock	○	●
Temperature monitor and homeothermic control module	FREE	FREE