


Part Number: **C30PU-RBD1443**



not to scale

Full Description: Catheter for rat bile duct, PU 3Fr 30cm, beveled tip, collars @ 0.7 and 5.0cm, sleeve. Fits 22ga.

Unit: each **Pouch Desc:** Rat BDC 3Fr PU 30cm **Pouch:** 4x10

Approx Dead Vol: 95 μ L **Connects to:** 22ga **EO Sterilized:**

Intended Species: rat **Vessel:** bile duct

Usage Notes: formerly PU-C30-BI01

TUBING

Tip: beveled **Length 1** **Length 2**
L1: 30 cm **L2:** cm
Material (IDxOD): PU 3Fr (.62x1.02 \pm .03mm) **Material (IDxOD):** cm

FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 0.7 cm	<input type="checkbox"/>
Feature B:	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 5.0 cm	<input type="checkbox"/>
Feature C:		c: cm	<input type="checkbox"/>
Proximal end			
Feature Z:	Silicone sleeve for 3Fr (2.2mm OD), 15mm long	z: 1.0 cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

PARTS ADDED TO POUCH

Part 1: **Part 2:**