


Part Number: **C10PU-MCA1301**



*not to scale*

**Full Description:** Catheter for mouse carotid artery, 1 to 3Fr, 8cm, collar @ 1.1cm. Fits 22ga.

**Unit:** each **Pouch Desc:** Mouse CAC PU 1-3Fr 8cm **Pouch:** 3x7.5

**Approx Dead Vol:** 17  $\mu$ L **Connects to:** 22ga **EO Sterilized:**

**Intended Species:** mouse **Vessel:** carotid artery

**Usage Notes:**

## TUBING

**Tip:** rounded **Length 1** **Length 2**  
**L1:** 2.5 cm **L2:** 5.5 cm  
**Material (IDxOD):** PU 1Fr (.18x0.35 $\pm$ .05mm) **Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm)

## FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 1Fr (0.6mm OD), 1mm wide	a: 1.1 cm	<input type="checkbox"/>
Feature B:		b: cm	<input type="checkbox"/>
Feature C:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	c: 2.5 cm	<input checked="" type="checkbox"/>
Proximal end			
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:**  **Part 2:**