


Part Number: **C30PU-RFV1418**



*not to scale*

**Full Description:** Catheter for rat femoral vein, PU 3Fr 30cm, collar @ 5.0cm. Fits 22ga.

**Unit:** each      **Pouch Desc:** Rat FVC 3Fr PU 30cm      **Pouch:** 4x10

**Approx Dead Vol:** 90  $\mu$ L      **Connects to:** 22ga      **EO Sterilized:**

**Intended Species:** rat      **Vessel:** femoral vein

**Usage Notes:**

## TUBING

**Tip:** rounded      *Length 1*      *Length 2*  
 L1: 30 cm      L2: cm  
**Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm)      **Material (IDxOD):**

## FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
<b>Feature A:</b>	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	a: 5.0 cm	<input type="checkbox"/>
<b>Feature B:</b>		b: cm	<input type="checkbox"/>
<b>Feature C:</b>		c: cm	<input type="checkbox"/>
<i>Proximal end</i>			
<b>Feature Z:</b>		z: cm	<input type="checkbox"/>
<b>Feature Y:</b>		y: cm	<input type="checkbox"/>
<b>Feature X:</b>		x: cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:**      **Part 2:**