

GEM- 800-125P Human Serum Albumin Powder

PRODUCT DESCRIPTION:

Human Albumin is prepared by cold-alcohol fractionation from pooled human plasma and lyophilized. MW ~66kDa.

Precautions: Non-Hazardous Ingredients: Blood Plasma (Human)

The product contains no substances that are considered hazardous to health at their given concentration. Although the human plasma used to manufacture this product has been tested and shown to be negative for infectious agents, no test is 100% accurate. All materials derived from human plasma should be handled as if capable of transmitting infection.

Not for clinical or therapeutic use.

INSTRUCTIONS FOR USE

Receipt and Storage:

Ships ambient. Store refrigerated at 2 to 8°C.

Appearance:

Cream Powder

Stability:

Expiration is lot dependent. If stored at 2-8°C, see Certificate of Analysis.

Microbiological Testing	
Microbial Growth	<2.99 cfu/mL
Viral Testing	
HBsAg	Non-Reactive
Anti-HCV	Non-Reactive
Anti-HIV-1/HIV-2	Non-Reactive
Syphilis	Negative
HIV-NAT	Not Detected
HCV-NAT	Not Detected
Physical Testing	
pH @ RT	Test and Report*
Purity	Test and Report
Endotoxin	≤2 EU/mg

*Based on a 5% Solution

For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.

QUALITY - TRACEABILITY - INTEGRITY

Life Science Production

The Science Centre, 1 Blaydon Road, Sandy, Bedfordshire, SG19 1RZ

Tel: +44 (0) 1234 889180 Email: sales@lifesciencegroup.co.uk Web: www.lifesciencegroup.co.uk

Life Science Production, Antibody Production Services and Diagnostic Kit Services are divisions of Life Science Group Ltd