

Biological Safety Cabinets Class II

NEOTERIC

Standard configuration:

- electric sockets
- three-section tabletop
- hinged front sash
- armrest
- pullout UV unit
- HEPA filters
- microprocessor control system
- pushbutton control panel
- frame stand



70% recirculation



Size (mm)	Dimensions (mm W x D x H)	Work Chamber (mm W x D x H)
1200	1200 x 770 x 2150	1105 x 610 x 700

SAVVY

Standard configuration:

- LED lighting
- electric sockets
- three-section tabletop
- hinged front sash
- armrest
- pullout UV unit
- HEPA filters
- microprocessor control system
- touchscreen control panel
- frame stand



70% recirculation




Size (mm)	Dimensions (mm W x D x H)	Work Chamber (mm W x D x H)
900	1000 x 770 x 2150	905 x 610 x 700
1200	1200 x 770 x 2150	1105 x 610 x 700
1500	1500 x 770 x 2150	1405 x 610 x 700
1800	1800 x 770 x 2150	1705 x 610 x 700

SAVVY SL

Standard configuration:

- LED lighting
- electric sockets
- three-section tabletop
- power-driven sliding front sash
- armrest
- pullout UV unit
- HEPA filters
- microprocessor control system
- touchscreen control panel
- frame stand
- voltage free output



70% recirculation



Size (mm)	Dimensions (mm W x D x H)	Work Chamber (mm W x D x H)
900	1000 x 770 x 2150	905 x 610 x 700
1200	1200 x 770 x 2150	1105 x 610 x 700
1500	1500 x 770 x 2150	1405 x 610 x 700
1800	1800 x 770 x 2150	1705 x 610 x 700

NEOTERIC B2 / SAVVY B2

Standard configuration:

- electric sockets
- three-section tabletop
- hinged front sash
- armrest
- pullout UV unit
- G4 filter
- HEPA filters
- microprocessor control system
- pushbutton/touchscreen control panel
- frame stand
- exhaust hood



70% recirculation



Size (mm)	Dimensions (mm W x D x H)	Work Chamber (mm W x D x H)
1200	1200 x 810 x 2325	1105 x 610 x 705
1500	1500 x 810 x 2220	1405 x 610 x 705
1800	1800 x 810 x 2325	1705 x 610 x 705

Biological Safety Cabinets Class III

PROTECT

Standard configuration:

- hinged front sash
- chamber gloves
- LED lighting
- removable electric sockets
- UV light in the work chamber
- G4 filter
- HEPA filters
- pass box
- microprocessor control system
- touchscreen control panel
- frame stand

Size (mm)	Dimensions (mm W x D x H)	Work Chamber (mm W x D x H)
1200	1720 x 750 x 1940	1140 x 640 x 675
1800	2320 x 750 x 1940	1740 x 640 x 675



PROTECT VIS-A-VIS

Standard configuration:

- hinged front sash
- chamber gloves
- LED lighting
- removable electric sockets
- UV light in the work chamber
- G4 filter
- HEPA filters
- pass box
- microprocessor control system
- touchscreen control panel
- frame stand

Size (mm)	Dimensions (mm W x D x H)	Work Chamber (mm W x D x H)
1200	1720 x 880 x 1940	1140 x 655 x 685
1800	2320 x 880 x 1940	1640 x 655 x 685



Vertical Laminar Flow Cabinets

LORICA

Standard configuration:

- electric sockets
- one-piece tabletop
- hinged front sash
- work opening screen
- UV light in the work chamber
- HEPA filters
- microprocessor control system
- pushbutton control panel
- frame stand



Size (mm)	Dimensions (mm W x D x H)	Work Chamber (mm W x D x H)
1200	1200 x 770 x 2150	1130 x 625 x 650
1500	1500 x 760 x 1870	1430 x 625 x 650
1800	1800 x 760 x 1870	1730 x 625 x 650

PCR Cabinets

PCR

Standard configuration:

- hinged front sash
- electric sockets
- one-piece tabletop
- pullout UV unit
- UV air recirculator
- frame stand



Size (mm)	Dimensions (mm W x D x H)	Work Chamber (mm W x D x H)
1000	1020 x 550 x 1317	1105 x 532 x 530