


Part Number: **C10PU-MFV1301**



*not to scale*

**Full Description:** Catheter for mouse femoral vein/artery, 1 to 3Fr, 10.5cm, collar @1.2cm. Fits 22ga.

**Unit:** each **Pouch Desc:** Mouse FVC PU 1-3Fr 10.5cm **Pouch:** 3x7.5

**Approx Dead Vol:** 16  $\mu$ L **Connects to:** 22ga **EO Sterilized:**

**Intended Species:** mouse **Vessel:** femoral vein or artery

**Usage Notes:**

## TUBING

**Tip:** rounded **Length 1** **Length 2**  
**L1:** 5 cm **L2:** 5 cm  
**Material (IDxOD):** PU 1Fr (.18x0.35 $\pm$ .05mm) **Material (IDxOD):** PU 3Fr (.62x1.02 $\pm$ .03mm)

## FEATURES

Distal (tip) end		Location	Fixed
Feature A:	Small silicone collar for 1Fr (0.6mm OD), 1mm wide	a: 1.2 cm	<input type="checkbox"/>
Feature B:		b: cm	<input type="checkbox"/>
Feature C:	Large silicone collar for 3Fr (1.7mm OD), 3mm wide	c: 5.0 cm	<input type="checkbox"/>
Proximal end			
Feature Z:		z: cm	<input type="checkbox"/>
Feature Y:		y: cm	<input type="checkbox"/>
Feature X:		x: cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:**  **Part 2:**